

Anti-MLKL Rabbit Polyclonal Antibody - Q8Nb16

Item Number: CM0023728



Introduction

MLKL Rabbit pAb is a polyclonal antibody targeting human MLKL (~54 kDa). It is applicable for WB, IF/ICC, IHC-P, and ELISA, with cross-reactivity to mouse and rat. Ideal for necroptosis and programmed cell death studies.

[Learn More](#)

Parameter	Value
Product Name	MLKL Rabbit pAb
Remarks/Alias	hMLKL; MLKL
Species	Human
Gene ID	197259
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 369-471 of human MLKL (NP_689862.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NB16
Protein Weight	54 kDa
Shipping	Ice bag