


Anti-MLL5/Kmt2E Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q8Izd2

Item Number: CM0028702



Introduction

Rabbit polyclonal antibody against human MLL5 (KMT2E). Validated for WB, IF-P, IHC-P, and ELISA. Reacts with human, mouse, and rat MLL5. Suitable for hematopoiesis, cell cycle, and chromatin regulation studies.

[Learn More](#)

Parameter	Value
Product Name	MLL5 Rabbit pAb
Remarks/Alias	MLL5; Nkp44L; ODLURO; SETD5B; HDCMC04P
Species	Human
GeneID	55904
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-500 of human MLL5 (NP_061152.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IZD2
Protein Weight	205kDa
Shipping	Ice bag