

Anti-MLLT1/ENL Polyclonal Antibody For Wb And Elisa - Q03111

Item Number: CM0022966



Introduction

Rabbit polyclonal antibody targeting human MLLT1/ENL, a YEATS domain chromatin reader. Validated for WB and ELISA; cross-reacts with human, mouse, and rat. Suitable for chromatin and leukemia research.

[Learn More](#)

Parameter	Value
Product Name	MLLT1 Rabbit pAb
Remarks/Alias	ENL; LTG19; YEATS1; MLLT1
Species	Human
Gene ID (Human)	4298
Gene ID	4298
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-270 of human MLLT1 (NP_005925.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03111
Protein Weight	62kDa
Shipping	Ice bag