

Anti-MLLT3/AF9 Rabbit Monoclonal Antibody For Wb, Elisa - P42568

Item Number: CM0005596



Introduction

High-quality Anti-MLLT3/AF9 Rabbit Monoclonal Antibody validated for WB and ELISA. Specific to human MLLT3/AF9, a chromatin reader in the super elongation complex essential for transcription regulation and hematopoietic stem cell maintenance.

[Learn More](#)

Parameter	Value
Product Name	MLLT3/AF9 Rabbit mAb
Remarks/Alias	AF9; YEATS3; MLLT3/AF9
Species	Human
GeneID (Human)	4300
GeneID	4300
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human MLLT3/AF9 (NP_004520.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P42568
Protein Weight	63kDa
Shipping	Ice bag