

Anti-MLLT3/AF9 Rabbit Polyclonal Antibody For Wb And Elisa - P42568

Item Number: CM0025649



Introduction

MLLT3/AF9 rabbit polyclonal antibody detects protein AF-9, a chromatin reader in the super elongation complex. Suitable for WB and ELISA in human samples, supporting epigenetics and leukemia research.

[Learn More](#)

Parameter	Value
Product Name	MLLT3/AF9 Rabbit pAb
Remarks/Alias	AF9; YEATS3; MLLT3/AF9
Species	Human
GeneID (Human)	4300
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human MLLT3/AF9 (NP_004520.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P42568
Protein Weight	63kDa
Shipping	Ice bag