


Anti-MLN Polyclonal Antibody For Wb, Ihc-P, Elisa - P12872

Item Number: CM0030170



Introduction

Anti-MLN polyclonal antibody raised in rabbit, validated for WB, IHC-P, ELISA. Recognizes human and mouse promotilin, a secreted hormone regulating gastrointestinal motility. Ideal for gastrointestinal research.

[Learn More](#)

Parameter	Value
Product Name	MLN Rabbit pAb
Remarks/Alias	MLN; motilin
Species	Human
GeneID (human)	4295
GeneID	4295
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-115 of human MLN (NP_002409.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P12872
Protein Weight	13kDa
Shipping	Ice bag