

Anti-Mlph Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bv36

Item Number: CM0025300

Introduction

Rabbit polyclonal antibody against human MLPH (Melanophilin), validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Useful for melanosome transport studies.

[Learn More](#)

Parameter	Value
Product Name	MLPH Rabbit pAb
Remarks/Alias	SLAC2-A; MLPH
Species	Human
Gene ID (Human)	79083
Gene ID	79083
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BV36
Protein Weight	66kDa
Shipping	Ice bag