


Anti-MLst8 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bvc4

Item Number: CM0016819



Introduction

KD-validated rabbit polyclonal antibody targeting MLST8 (GbetaL, mLST8), suitable for Western blot and ELISA in human, mouse, and rat. Recognizes the mTOR complex subunit MLST8, essential for cell growth and survival signaling.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] MLST8 Rabbit pAb
Remarks/Alias	GBL; LST8; POP3; WAT1; GbetaL; [KD Validated] MLST8
Species	Human
GeneID (Human)	64223
GeneID	64223
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-326 of human MLST8 (NP_001186103.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BVC4
Protein Weight	36kDa
Shipping	Ice bag