


Anti-Mlycd Polyclonal Antibody For Wb, Elisa - 095822

Item Number: CM0021201



Introduction

Rabbit polyclonal antibody targeting human MLYCD (Malonyl-CoA decarboxylase), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for fatty acid metabolism and mitochondrial research.

[Learn More](#)

Parameter	Value
Product Name	MLYCD Rabbit pAb
Remarks/Alias	MCD; MLYCD
Species	Human
GeneID (Human)	23417
GeneID	23417
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-190 of human MLYCD (NP_036345.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	095822
Protein Weight	55kDa
Shipping	Ice bag