


Anti-Mmachc Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y4U1

Item Number: CM0017560



Introduction

Rabbit polyclonal antibody against human MMACHC (cbIC), validated for Western blot and ELISA. Cross-reacts with mouse and rat. For cobalamin metabolism and methylation pathway studies.

[Learn More](#)

Parameter	Value
Product Name	MMACHC Rabbit pAb
Remarks/Alias	cbIC; MMACHC
Species	Human
GeneID (Human)	25974
GeneID	25974
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 63-282 of human MMACHC (NP_056321.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4U1
Protein Weight	32kDa
Shipping	Ice bag