


Anti-Mmp1 Rabbit Monoclonal Antibody For Wb, Elisa - P03956

Item Number: CM0006753



Introduction

Recombinant rabbit monoclonal antibody targeting human MMP1 (collagenase), suitable for WB and ELISA. Cross-reacts with human, mouse, rat. Key reagent for extracellular matrix degradation research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	MMP1 Rabbit mAb
Remarks/Alias	CLG; CLGN; MMP1
Species	Human
GeneID (Human)	4312
GeneID	4312
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 93-210 of human MMP1 (NP_002412.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P03956
Protein Weight	54kDa
Shipping	Ice bag