


Anti-Mmp10 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09238

Item Number: CM0026594



Introduction

Rabbit polyclonal antibody against human MMP10 (Stromelysin-2, SL-2, STMY2). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	MMP10 Rabbit pAb
Remarks/Alias	SL-2; STMY2; MMP10
Species	Human
GeneID (human)	4319
GeneID	4319
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 237-476 of human MMP10 (NP_002416.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09238
Protein Weight	54kDa
Shipping	Ice bag