

Anti-Mmp11 Rabbit Polyclonal Antibody For Wb, Elisa - P24347

Item Number: CM0028467

Introduction

Rabbit polyclonal antibody against human MMP11 (stromelysin-3, matrix metalloproteinase-11) for Western blot and ELISA. Reacts with human, mouse, rat. Targets amino acids 239-488. Implicated in epithelial malignancy progression.

[Learn More](#)

Parameter	Value
Product Name	MMP11 Rabbit pAb
Remarks/Alias	ST3; SL-3; STMY3; MMP11
Species	Human
GeneID (human)	4320
GeneID	4320
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 239-488 of human MMP11 (NP_005931.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P24347
Protein Weight	55kDa
Shipping	Ice bag