


Anti-Mmp12 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P39900

Item Number: CM0023831



Introduction

MMP12 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA in human, mouse, rat. Targets macrophage metalloelastase, involved in tissue remodeling and elastolysis. Ideal for inflammation and ECM studies.

[Learn More](#)

Parameter	Value
Product Name	MMP12 Rabbit pAb
Remarks/Alias	ME; HME; MME; MMP-12; MMP12
Species	Human
Gene ID (Human)	4321
Gene ID	4321
Immunogen	Synthetic peptide (amino acids 200-300 of human MMP12, NP_002417.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P39900
Protein Weight	54kDa
Shipping	Ice bag