


Anti-Mmp13 Monoclonal Antibody For Wb, Elisa - P45452

Item Number: CM0008001



Introduction

Rabbit monoclonal MMP13 antibody validated for WB, ELISA. Recognizes human MMP13 (collagenase 3), SWISS P45452, with cross-reactivity to mouse. Suitable for extracellular matrix degradation and bone development studies.

[Learn More](#)

Parameter	Value
Product Name	MMP13 Rabbit mAb
Remarks/Alias	CLG3; MDST; MANDP1; MMP-13
Species	Human
GeneID (human)	4322
GeneID	4322
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 104-272 of human MMP13 (NP_002418.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P45452
Protein Weight	54 kDa
Shipping	Ice bag