


Anti-Mmp13 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P45452

Item Number: CM0002750



Introduction

High-quality rabbit monoclonal antibody targeting human MMP13 (collagenase 3). Validated for WB, IF/ICC, and ELISA. Suitable for extracellular matrix degradation, bone development, and cancer invasion studies.

[Learn More](#)

Parameter	Value
Product Name	MMP13 Rabbit mAb
Remarks/Alias	CLG3; MDST; MANDP1; MMP-13; MMP13
Species	Human
GeneID	4322
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 350-450 of human MMP13 (P45452))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P45452
Protein Weight	54kDa
Shipping	Ice bag