


Anti-Mmp14 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - P50281

Item Number: CM0025669



Introduction

MMP14 Rabbit Polyclonal Antibody validated for WB, IF/ICC, IF-P, ELISA in human and mouse. Detects a 66kDa membrane-type matrix metalloproteinase involved in ECM remodeling and cell migration.

[Learn More](#)

Parameter	Value
Product Name	MMP14 Rabbit pAb
Remarks/Alias	MMP-14; MMP-X1; MT-MMP; MT1MMP; MTMMP1; WNCHRS; MT1-MMP; MT-MMP 1
Species	Human
GeneID (human)	4323
GeneID	4323
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 21-538 of human MMP14 (NP_004986.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P50281
Protein Weight	66kDa
Shipping	Ice bag