


Anti-Mmp14/Mt1-Mmp Monoclonal Antibody For Wb, Ihc-P, Fc, Elisa - P50281

Item Number: CM0003304



Introduction

High-quality rabbit monoclonal antibody targeting human MMP14/MT1-MMP, validated for WB, IHC-P, FC, and ELISA. Cross-reacts with mouse and rat. Ideal for extracellular matrix and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MMP14/MT1-MMP Rabbit mAb
Remarks/Alias	MMP-14; MMP-X1; MT-MMP; MT1MMP; MTMMP1; WNCHRS; MT1-MMP; MT-MMP 1; MMP14/MT1-MMP
Species	Human
Gene ID (Human)	4323
Gene ID	4323
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human MMP14 (P50281)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, FC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50281
Protein Weight	66kDa
Shipping	Ice bag