


Anti-Mmp14/Mt1-Mmp Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P50281

Item Number: CM0030110



Introduction

Rabbit polyclonal antibody targeting MMP14/MT1-MMP, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat MMP14, suitable for investigating collagenolysis, cell migration, and Notch signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	MMP14/MT1-MMP Rabbit pAb
Remarks/Alias	MMP-14; MMP-X1; MT-MMP; MT1MMP; MTMMP1; WNCHRS; MT1-MMP; MT-MMP 1; MMP14/MT1-MMP
Species	Human
GeneID (Human)	4323
GeneID	4323
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50281
Protein Weight	66kDa
Shipping	Ice bag