


Anti-Mmp16 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P51512

Item Number: CM0016509



Introduction

High-quality rabbit polyclonal antibody against human MMP16 (MT3-MMP), suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse. Target is a membrane-type matrix metalloproteinase involved in extracellular matrix degradation and vascular remodeling.

[Learn More](#)

Parameter	Value
Product Name	MMP16 Rabbit pAb
Remarks/Alias	MMP-X2; C8orf57; MT-MMP2; MT-MMP3; MT3-MMP; MMP16
Species	Human
GeneID (Human)	4325
GeneID	4325
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 290-565 of human MMP16 (NP_005932.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51512
Protein Weight	70kDa
Shipping	Ice bag