


Anti-Mmp2 Rabbit Monoclonal Antibody For Wb, Elisa - P08253

Item Number: CM0002812



Introduction

Anti-MMP2 rabbit monoclonal antibody validated for Western blot and ELISA, reacting with human, mouse, and rat MMP-2, a secreted metalloproteinase involved in angiogenesis and tumor invasion.

[Learn More](#)

Parameter	Value
Product Name	MMP2 Rabbit mAb
Remarks/Alias	CLG4; MONA; CLG4A; MMP-2; TBE-1; MMP-II; MMP2
Species	Human
GeneID (Human)	4313
GeneID	4313
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 500-600 of human MMP2 (P08253).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08253
Protein Weight	74kDa
Shipping	Ice bag