


Anti-Mmp2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08253

Item Number: CM0014797



Introduction

This MMP2 rabbit polyclonal antibody detects human, mouse, rat MMP2, validated for WB, IF/ICC, ELISA. Suitable for studying angiogenesis, tumor invasion, and tissue remodeling.

[Learn More](#)

Parameter	Value
Product Name	MMP2 Rabbit pAb
Remarks/Alias	CLG4; MONA; CLG4A; MMP-2; TBE-1; MMP-II; MMP2
Species	Human
GeneID (Human)	4313
Gene ID	4313
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 40-280 of human MMP2 (NP_004521.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08253
Protein Weight	74kDa
Shipping	Ice bag