

Anti-Mmp20 Rabbit Polyclonal Antibody For Wb, Elisa - O60882

Item Number: CM0029160



Introduction

MMP20 Rabbit Polyclonal Antibody validated for WB and ELISA, reacts with Human, Mouse, Rat. Target is enamel metalloproteinase involved in tooth enamel formation. Ideal for dental and skeletal research.

[Learn More](#)

Parameter	Value
Product Name	MMP20 Rabbit pAb
Remarks/Alias	AI2A2; MMP-20; MMP20
Species	Human
Gene ID (Human)	9313
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 108-310 of human MMP20 (NP_004762.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60882
Protein Weight	54kDa
Shipping	Ice bag