

Anti-Mmp25 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Npa2

Item Number: CM0022553

Introduction

Rabbit polyclonal antibody against human MMP25. Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for protease and extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	MMP25 Rabbit pAb
Remarks/Alias	MMP20; MMPL1; MMP-25; MMP20A; MT6MMP; MTMMP6; MT-MMP6; MT6-MMP; MT-MMP 6; MMP25
Species	Human
Gene ID	64386
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 321-539 of human MMP25 (NP_071913.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPA2
Protein Weight	63kDa
Shipping	Ice bag