

Anti-Mmp25 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Npa2

Item Number: CM0027427

Introduction

Rabbit polyclonal antibody targeting human MMP25 (MT6-MMP/leukolysin). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Recombinant immunogen (aa 350-540). ~63 kDa. For cancer and leukocyte research.

[Learn More](#)

Parameter	Value
Product Name	MMP25 Rabbit pAb
Remarks/Alias	MMP20; MMPL1; MMP-25; MMP20A; MT6MMP; MTMMP6; MT-MMP6; MT6-MMP; MT-MMP 6; MMP25
Species	Human
GeneID (Human)	64386
GeneID	64386
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-540 of human MMP25 (NP_071913.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPA2
Protein Weight	63kDa
Shipping	Ice bag