

Anti-Mmp3 Polyclonal Antibody For Wb, Elisa - P28862

Item Number: CM0028979

Introduction

Rabbit polyclonal antibody against MMP3 (Stromelysin-1) for WB and ELISA applications. Recognizes mouse MMP3 with a protein weight of 54kDa. Suitable for extracellular matrix degradation and innate immunity research.

[Learn More](#)

Product Name	MMP3 Rabbit pAb
Remarks/Alias	SL-1; STMY; STR1; CHDS6; MMP-3; STMY1
Species	Mouse
GeneID (Human)	4314
GeneID	17392
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 385-477 of human MMP3 (NP_034939.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P28862
Protein Weight	54kDa
Shipping	Ice bag