


Anti-Mmp3 (Stromelysin-1) Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08254

Item Number: CM0006706



Introduction

Rabbit monoclonal antibody targeting human MMP3 (Stromelysin-1). Validated for WB, IHC-P, IF/ICC, and ELISA with cross-reactivity to mouse and rat. Ideal for ECM remodeling and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	MMP3 Rabbit mAb
Remarks/Alias	SL-1; STMY; STR1; CHDS6; MMP-3; STMY1; MMP3
Species	Human
Gene ID	4314 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 385-477 of human MMP3 (NP_002413.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08254
Protein Weight	54kDa
Shipping	Ice bag