


Anti-Mmp3 (Stromelysin-1) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P08254

Item Number: CM0003788



Introduction

Rabbit monoclonal antibody against human MMP3 (Stromelysin-1). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with rat. Ideal for ECM degradation and inflammatory disease research.

[Learn More](#)

Parameter	Value
Product Name	MMP3 Rabbit mAb
Remarks/Alias	SL-1; STMY; STR1; CHDS6; MMP-3; STMY1; MMP3
Species	Human
Gene ID (Human)	4314
Gene ID	4314
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 385-477 of human MMP3 (NP_002413.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P08254
Protein Weight	54 kDa
Shipping	Ice bag