


Anti-Mmp7 (Matrilysin) Rabbit Monoclonal Antibody For Wb, Elisa - P09237

Item Number: CM0003948



Introduction

MMP7 rabbit monoclonal antibody for WB and ELISA. Specifically reacts with human matrilysin (MMP-7), ~30 kDa. Ideal for extracellular matrix remodeling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	MMP7 Rabbit mAb
Remarks/Alias	MMP-7; MPSL1; PUMP-1; MMP7
Species	Human
Gene ID (Human)	4316
Gene ID	4316
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-267 of human MMP7 (P09237)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P09237
Protein Weight	30kDa
Shipping	Ice bag