


Anti-Mmp7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P09237

Item Number: CM0017911



Introduction

High-quality rabbit polyclonal antibody targeting human MMP7 (Matrilysin/PUMP-1), validated for Western blotting, IHC-P, IF/ICC, and ELISA. Cross-reacts with rat. Ideal for studying matrix metalloproteinase-7 in cancer, tissue remodeling, and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	MMP7 Rabbit pAb
Remarks/Alias	MMP-7; MPSL1; PUMP-1; MMP7
Species	Human
GeneID (Human)	4316
GeneID	4316
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human MMP7 (NP_002414.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	P09237
Protein Weight	30kDa
Shipping	Ice bag