


Anti-Mmp8 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P22894

Item Number: CM0020434



Introduction

Rabbit polyclonal antibody targeting MMP8 (neutrophil collagenase) for WB, IHC-P, ELISA applications. Reacts with human, mouse, rat. Suitable for extracellular matrix and inflammation research. SWISS: P22894.

[Learn More](#)

Parameter	Value
Product Name	MMP8 Rabbit pAb
Remarks/Alias	HNC; CLG1; MMP-8; PMNL-CL; MMP8
Species	Human
Gene ID (Human)	4317
Gene ID	4317
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 101-467 of human MMP8 (NP_002415.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22894
Protein Weight	53kDa
Shipping	Ice bag