

Anti-Mmp9 Polyclonal Antibody For Wb, Elisa - P41245

Item Number: CM0020291

Introduction

High-quality rabbit polyclonal antibody against mouse MMP9 (gelatinase B, 92 kDa type IV collagenase). Validated for WB and ELISA, ideal for investigating ECM remodeling, bone resorption, and metastasis in mouse research.

[Learn More](#)

Parameter	Value
Product Name	MMP9 Rabbit pAb
Remarks/Alias	Clg4b, Gel B, MMP-9, B/MMP9, pro-MMP-9, Pro-MMP9
Species	Mouse
GeneID (Human)	4318
GeneID (Mouse)	17395
Immunogen	Recombinant protein corresponding to a sequence within amino acids 108-730 of mouse MMP9 (NP_038627.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P41245
Protein Weight	80KDa
Shipping	Ice bag