


Anti-Mmp9 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P14780

Item Number: CM0016196



Introduction

Rabbit polyclonal antibody against human MMP9 (gelatinase B, CLG4B). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. For extracellular matrix and leukocyte migration research.

[Learn More](#)

Parameter	Specification
Product Name	MMP9 Rabbit pAb
Remarks/Alias	GELB; CLG4B; MMP-9; MANDP2; MMP9
Species	Human
GeneID (Human)	4318
GeneID	4318
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 538-707 of human MMP9 (NP_004985.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14780
Protein Weight	78 kDa
Shipping	Ice bag