

Anti-Mmut Monoclonal Antibody For Wb, Ihc-P, Elisa - P22033

Item Number: CM0012352



Introduction

High-quality anti-MMUT rabbit monoclonal antibody validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for studying methylmalonyl-CoA mutase function in metabolic research.

[Learn More](#)

Parameter	Value
Product Name	MMUT Rabbit mAb
Remarks/Alias	MCM; MUT; MMUT
Species	Human
GeneID	4594
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human MMUT (P22033).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22033
Protein Weight	83kDa
Shipping	Ice bag