

Anti-Mnat1 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P51948

Item Number: CM0019080



Introduction

Rabbit polyclonal antibody against human MNAT1 (MAT1/CDK7 assembly factor), suitable for WB, IP, ELISA. Cross-reacts with human and mouse. For cell cycle and transcription regulation research. Protein weight 36 kDa.

[Learn More](#)

Parameter	Value
Product Name	MNAT1 Rabbit pAb
Remarks/Alias	MAT1; TFB3; CAP35; RNF66; MNAT1
Species	Human
Gene ID	4331
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-309 of human MNAT1 (NP_002422.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51948
Protein Weight	36kDa
Shipping	Ice bag