


Anti-Mob1A/B Monoclonal Antibody For Wb, Elisa - Q9H8S9 / Q7L9L4

Item Number: CM0004902



Introduction

High-quality rabbit monoclonal antibody targeting human MOB1A/B, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for Hippo signaling and cancer research. One-stop sourcing from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	MOB1A/B Rabbit mAb
Remarks/Alias	MOB1; MATS1; Mob4B; C2orf6; MOBK1B; MABKL1B; MOBKL1B; MOB1A/B
Species	Human
GeneID (Human)	55233/92597
GeneID	55233/92597
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human MOB1A/B (NP_060691.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H8S9 / Q7I9L4
Protein Weight	25 kDa
Shipping	Ice bag