


Anti-Mob1A/Mob1B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H8S9 / Q7L9L4

Item Number: CM0021390



Introduction

Rabbit polyclonal antibody targeting MOB1A/MOB1B, validated for WB, IF/ICC, ELISA with human, mouse, rat reactivity. MOB1 proteins are Hippo pathway activators regulating organ size and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	MOB1A/MOB1B Rabbit pAb
Remarks/Alias	MOB1A/MOB1B
Species	Human
Gene ID (Human)	55233/92597
Gene ID	55233/92597
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H8S9/Q7I9L4
Shipping	Ice bag