


Anti-Mob4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y3A3

Item Number: CM0024968



Introduction

Rabbit polyclonal antibody targeting human MOB4, validated for Western blot, IHC-paraffin, and ELISA. Cross-reacts with mouse and rat. Useful for investigating STRIPAK complex functions in cell signaling and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	MOB4 Rabbit pAb
Remarks/Alias	2C4D; MOB1; MOB3; PHOCN; PREI3; CGI-95; MOBKL3; MOB4
Species	Human
Gene ID (Human)	25843
Gene ID	25843
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human MOB4 (NP_056202.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3A3
Protein Weight	26kDa
Shipping	Ice bag