

Anti-Mobp Rabbit Polyclonal Antibody For Wb And Elisa - Q13875

Item Number: CM0020795

Introduction

High-quality anti-MOBP rabbit polyclonal antibody validated for WB and ELISA applications. Reacts with human, mouse, and rat MOBP, a myelin-associated protein involved in myelin sheath compaction. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	MOBP Rabbit pAb
Remarks/Alias	MOBP
Species	Human
Gene ID (Human)	4336
Gene ID	4336
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human MOBP (NP_891980.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13875
Protein Weight	21kDa
Shipping	Ice bag