

Anti-Moesin Monoclonal Antibody For Wb, Ihc-P, Elisa - P26038

Item Number: CM0009165



Introduction

Rabbit monoclonal antibody against human Moesin, validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for studying ERM family proteins, cell cortex organization, and immune signaling. UniProt P26038.

[Learn More](#)

Parameter	Value
Product Name	Moesin Rabbit mAb
Remarks/Alias	HEL70; IMD50; Moesin
Species	Human
Gene ID (Human)	4478
Gene ID	4478
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 400-500 of human MSN (P26038).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26038
Protein Weight	68kDa
Shipping	Ice bag