

Anti-Moesin Polyclonal Antibody For Wb, Elisa - P26038

Item Number: CM0024417



Introduction

Rabbit polyclonal antibody against human moesin (MSN), validated for WB and ELISA. Cross-reacts with human and mouse. Ideal for studying ERM family proteins and actin cytoskeleton regulation.

[Learn More](#)

Parameter	Value
Product Name	Moesin Rabbit pAb
Remarks/Alias	HEL70; IMD50; Moesin
Species	Human
GeneID (Human)	4478
GeneID	4478
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human Moesin (NP_002435.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P26038
Protein Weight	68 kDa
Shipping	Ice bag