

Anti-Moesin Polyclonal Antibody For Wb, If/Icc, Elisa - P26038

Item Number: CM0015088



Introduction

Rabbit polyclonal antibody targeting human Moesin (MSN), validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Ideal for studying ERM protein functions in cell cortex organization, immune regulation, and cytoskeletal dynamics.

[Learn More](#)

Parameter	Value
Product Name	Moesin Rabbit pAb
Remarks/Alias	HEL70; IMD50; Moesin
Species	Human
GeneID (Human)	4478
GeneID	4478
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 458-557 of human Moesin (NP_002435.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26038
Protein Weight	68kDa
Shipping	Ice bag