

## Anti-Mogs Rabbit Polyclonal Antibody For Wb And Elisa - Q13724

Item Number: CM0013273

### Introduction

This rabbit polyclonal antibody detects human MOGS (mannosyl-oligosaccharide glucosidase, DER7/GCS1), ~92kDa, for WB and ELISA. Cross-reacts with mouse and rat. Excellent for N-glycan degradation research - Q13724.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | MOGS Rabbit pAb                                                                                                   |
| Remarks/Alias    | DER7; GCS1; CDG2B; CWH41; MOGS                                                                                    |
| Species          | Human                                                                                                             |
| Gene ID (Human)  | 7841                                                                                                              |
| Gene ID          | 7841                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 60-320 of human MOGS (NP_006293.2). |
| Source           | Rabbit                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                             |
| Application      | WB, ELISA                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                 |
| SWISS            | Q13724                                                                                                            |
| Protein Weight   | 92kDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |