

Anti-Mok Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uq07

Item Number: CM0022938

Introduction

Rabbit polyclonal antibody against human MOK (RAGE-1), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for kinase signaling and ciliogenesis research.

[Learn More](#)

Parameter	Value
Product Name	MOK Rabbit pAb
Remarks/Alias	RAGE; RAGE1; STK30; RAGE-1; MOK
Species	Human
Gene ID (Human)	5891
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human MOK (NP_055041.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQ07
Protein Weight	48kDa
Shipping	Ice bag