

Anti-Monoamine Oxidase A (Maoa) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P21397

Item Number: CM0009662



Introduction

Rabbit monoclonal antibody against human Monoamine Oxidase A (MAOA), validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat. Ideal for neurotransmitter metabolism and mitochondrial function research.

[Learn More](#)

Parameter	Value
Product Name	Monoamine Oxidase A (MAOA) Rabbit mAb
Remarks/Alias	BRNRS; MAO-A; Monoamine Oxidase A (MAOA)
Species	Human
GeneID (Human)	4128
GeneID	4128
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 428-527 of human Monoamine Oxidase A (MAOA) (P21397).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21397
Protein Weight	60kDa
Shipping	Ice bag