

Anti-Monoamine Oxidase A (Maoa) Rabbit Polyclonal Antibody For If/Icc, Elisa - P21397

Item Number: CM0017498



Introduction

High-quality rabbit polyclonal antibody targeting human Monoamine Oxidase A (MAOA), validated for IF/ICC and ELISA. Detects mitochondrial MAOA, a key enzyme in neurotransmitter degradation. Suitable for neuroscience and cardiovascular research.

[Learn More](#)

Parameter	Value
Product Name	Monoamine Oxidase A (MAOA) Rabbit pAb
Remarks/Alias	BRNRS; MAO-A; Monoamine Oxidase A (MAOA)
Species	Human
GeneID (human)	4128
GeneID	4128
Immunogen	Recombinant fusion protein containing amino acids 1-260 of human Monoamine Oxidase A (NP_000231.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P21397
Protein Weight	61kDa
Shipping	Ice bag