

## Anti-Monoamine Oxidase A (Maoa) Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P21397

Item Number: CM0015490



### Introduction

Rabbit polyclonal antibody against human monoamine oxidase A (MAOA) for Western blot, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for neurotransmitter metabolic studies.

[Learn More](#)

| Parameter        | Value                                                                                                                          |
|------------------|--------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Monoamine Oxidase A (MAOA) Rabbit pAb                                                                                          |
| Remarks/Alias    | BRNRS; MAO-A; Monoamine Oxidase A (MAOA)                                                                                       |
| Species          | Human                                                                                                                          |
| GeneID (human)   | 4128                                                                                                                           |
| GeneID           | 4128                                                                                                                           |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human Monoamine Oxidase A (NP_000231.1) |
| Source           | Rabbit                                                                                                                         |
| Category         | Polyclonal Antibodies                                                                                                          |
| Application      | WB, IF/ICC, IHC-P, ELISA                                                                                                       |
| Cross Reactivity | Human, Mouse, Rat                                                                                                              |
| SWISS            | P21397                                                                                                                         |
| Protein Weight   | 61 kDa                                                                                                                         |
| Shipping         | Ice bag                                                                                                                        |