


Anti-Monoamine Oxidase B (Maob) Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P27338

Item Number: CM0003799



Introduction

High-quality rabbit monoclonal antibody targeting human Monoamine Oxidase B (MAOB). Validated for WB, IHC-P, ELISA with cross-reactivity to mouse, rat. Suitable for studying neurotransmitter metabolism and neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	Monoamine Oxidase B (MAOB) Rabbit mAb
Remarks/Alias	MAOB; monoamine oxidase B; Monoamine Oxidase B (MAOB)
Species	Human
GeneID (Human)	4129
GeneID	4129
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Monoamine Oxidase B (MAOB) (P27338).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27338
Protein Weight	59kDa
Shipping	Ice bag