

Anti-Morf4L2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q15014

Item Number: CM0028878

Introduction

Rabbit polyclonal antibody targeting human MORF4L2, validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying NuA4 complex, chromatin regulation, and DNA repair pathways.

[Learn More](#)

Parameter	Value
Product Name	MORF4L2 Rabbit pAb
Remarks/Alias	MRGX, MORFL2, MORF4L2
Species	Human
Gene ID (Human)	9643
Gene ID	9643
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-282 of human MORF4L2 (NP_001135890.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15014
Protein Weight	32kDa
Shipping	Ice bag