

Abflo 647 Anti-Mouse Cd107B (Lamp2) Monoclonal Antibody For Flow Cytometry - P17047

Item Number: CM0008198



Introduction

ABflo 647-conjugated rabbit monoclonal antibody targeting mouse CD107b (LAMP2) validated for flow cytometry. Reacts with mouse. Ideal for investigating lysosome biogenesis, autophagy, and antigen presentation pathways.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD107b mAb
Remarks/Alias	Mac3; LGP-B; CD107b; Lamp-2; Lamp II; Lamp-2a; Lamp-2b; Lamp-2c
Species	Mouse
GeneID (Human)	3920
GeneID	16784
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-189 of mouse CD107b (NP_001017959.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P17047
Protein Weight	46kDa
Shipping	Ice bag